

BRNO UNIVERSITY OF TECHNOLOGY

Faculty of Information Technology

Detection, Extraction and Measurement of the Contour and Circumference of the Metacarpal
Bones in X-rays of the Human Hand

Matej Otčenáš

1 Namespace Index	1
1.1 Packages	1
2 Class Index	2
2.1 Class List	2
3 Namespace Documentation	2
3.1 main Namespace Reference	2
3.1.1 Detailed Description	2
3.1.2 Function Documentation	2
3.2 Modules.canny Namespace Reference	2
3.2.1 Detailed Description	3
3.2.2 Function Documentation	3
3.3 Modules.config Namespace Reference	3
3.3.1 Detailed Description	3
3.3.2 Function Documentation	3
3.3.3 Variable Documentation	5
3.4 Modules.roi Namespace Reference	5
3.4.1 Detailed Description	5
3.4.2 Function Documentation	5
3.5 Modules.train Namespace Reference	5
3.5.1 Detailed Description	6
3.6 run Namespace Reference	6
3.6.1 Detailed Description	6
3.6.2 Variable Documentation	6
3.7 setup Namespace Reference	7
3.7.1 Detailed Description	7
3.7.2 Variable Documentation	7
4 Class Documentation	7
4.1 Modules.train.Trainer Class Reference	7
4.1.1 Detailed Description	8
4.1.2 Constructor & Destructor Documentation	8
4.1.3 Member Function Documentation	9
4.1.4 Member Data Documentation	9
Index	11

1 Namespace Index

1.1 Packages

Here are the packages with brief descriptions (if available):

main	2
Modules.canny	2
Modules.config	3
Modules.roi	5
Modules.train	5
run	6
setup	7

2 Class Index

2.1 Class List

Here are the classes, structs, unions and interfaces with brief descriptions:

Modules.train.Trainer	7
---------------------------------------	---

3 Namespace Documentation

3.1 main Namespace Reference

Functions

- def [main](#) ()

3.1.1 Detailed Description

`main.py`: Main module for running the contour detection and bone measurement algorithm.

3.1.2 Function Documentation

3.1.2.1 `main()` `def main.main ()`

Main function starts an algorithm

3.2 Modules.canny Namespace Reference

Functions

- def [Canny](#) (img_mask)

3.2.1 Detailed Description

canny.py: Module uses basic approach for contour extraction using OpenCV functions such as threshold or Canny edge detection algorithm.

3.2.2 Function Documentation

3.2.2.1 Canny() `def Modules.canny.Canny (`
 img_mask)

Method computes edges of given image mask using right threshold and Canny edge detection algorithm.

Parameters

img_mask : list
 Input image mask

Returns

array
 matrix of an image

3.3 Modules.config Namespace Reference

Functions

- `def TestMeta ()`
- `def TestROI ()`
- `def TrainMeta ()`
- `def TrainROI ()`

Variables

- `int MODEL_USE1 = 1`
- `int MODEL_USE2 = 1`

3.3.1 Detailed Description

config.py: Configuration module for training and inference. Values are specific for use of Detectron2 library.
https://detectron2.readthedocs.io/en/latest/tutorials/getting_started.html
<https://github.com/facebookresearch/detectron2>

NOTE:

Models were primary trained on Google Colab(<https://colab.research.google.com/>) due to high quality of GPU utilization.

3.3.2 Function Documentation

3.3.2.1 TestMeta() `def Modules.config.TestMeta ()`

Function for inference on given dataset of hands based on the trained model. Coming from 'TrainMeta()' function.

Returns

str, str, str, str, int, int, float
set of multiple parameters for model prediction

3.3.2.2 TestROI() `def Modules.config.TestROI ()`

Function for inference on given dataset of ROI of third metacarpal bone based on the trained model. Coming from 'TrainROI()' function.

Returns

str, str, str, str, int, int, float
set of multiple parameters for model prediction

3.3.2.3 TrainMeta() `def Modules.config.TrainMeta ()`

Function for training custom model based on given image annotations using pretrained deep neural network (DNN) called Mask RCNN for instance segmentation. Models utilize COCO(common objects in context) large-scale detection dataset. Training is provided on full x-ray image of human hand, where the third metacarpal is the only important.

Returns

str, str, str, str, int, int, int, float, int, int
set of multiple parameters for model configuration

3.3.2.4 TrainROI() `def Modules.config.TrainROI ()`

Function for training custom model based on given image annotations using pretrained deep neural network (DNN) called Mask RCNN for instance segmentation. Models utilize COCO(common objects in context) large-scale detection dataset. Training is provided on region of interest (ROI) of detected bones where the bone width is shortest.

Returns

str, str, str, str, int, int, int, float, int, int
set of multiple paremeters for model configuration

3.3.3 Variable Documentation

3.3.3.1 MODEL_USE1 `int Modules.config.MODEL_USE1 = 1`

3.3.3.2 MODEL_USE2 `int Modules.config.MODEL_USE2 = 1`

3.4 Modules.roi Namespace Reference

Functions

- def [ROI](#) (boxes)

3.4.1 Detailed Description

`roi.py`: Module obtains bounding box points coordinates.

3.4.2 Function Documentation

3.4.2.1 ROI() `def Modules.roi.ROI (
 boxes)`

Function extracts top left and bottom right point from given bounding box.

Parameters

boxes: array

Bounding box of processed bone.

Returns

int, int, int, int

set of coordinates as separate numbers

3.5 Modules.train Namespace Reference

Classes

- class [Trainer](#)

3.5.1 Detailed Description

`train.py`: Module for training the model based on custom dataset (annotated metacarpal bones or annotated ROI).

3.6 run Namespace Reference

Variables

- `ap` = `argparse.ArgumentParser()`
- string `arg` = `'python -W ignore main.py'`
- `args` = `vars(ap.parse_args())`
- `help`
- `input_dir` = `args["input"]`
- `output_dir` = `args["output"]`
- `required`
- `start_image` = `args["name"]`

3.6.1 Detailed Description

`run.py`: This module starts the entire third metacarpal edge detection and measurement program.

3.6.2 Variable Documentation

3.6.2.1 `ap` `run.ap = argparse.ArgumentParser()`

3.6.2.2 `arg` string `run.arg = 'python -W ignore main.py'`

3.6.2.3 `args` `run.args = vars(ap.parse_args())`

3.6.2.4 `help` `run.help`

3.6.2.5 `input_dir` `run.input_dir = args["input"]`

3.6.2.6 output_dir `run.output_dir = args["output"]`

3.6.2.7 required `run.required`

3.6.2.8 start_image `run.start_image = args["name"]`

3.7 setup Namespace Reference

Variables

- `ext_modules`
- `name`

3.7.1 Detailed Description

`setup.py`: Setup module for creating shared object libraries and C files for precompiling the program to faster

3.7.2 Variable Documentation

3.7.2.1 ext_modules `setup.ext_modules`

3.7.2.2 name `setup.name`

4 Class Documentation

4.1 Modules.train.Trainer Class Reference

Public Member Functions

- `def __init__` (self, `COCO_NAME`, `COCO_ANNOTS`, `MODEL_TYPE`, `WEIGHT_PATH`, `MAX_ITER_META`, `NUM_WORKERS`, `IMS_PER_BATCH`, `BASE_LR`, `BATCH_SIZE_PER_IMAGE`, `NUM_CLASSES`)
- `def run` (self)

Public Attributes

- [BASE_LR](#)
- [BATCH_SIZE_PER_IMAGE](#)
- [COCO_ANNOTS](#)
- [COCO_NAME](#)
- [IMS_PER_BATCH](#)
- [MAX_ITER_META](#)
- [MODEL_TYPE](#)
- [NUM_CLASSES](#)
- [NUM_WORKERS](#)
- [WEIGHT_PATH](#)

4.1.1 Detailed Description

Description

Trainer class including method for training the model.

Methods

`run()`
 Starts the training process

4.1.2 Constructor & Destructor Documentation

4.1.2.1 `__init__()` `def Modules.train.Trainer.__init__ (`
 `self,`
 `COCO_NAME,`
 `COCO_ANNOTS,`
 `MODEL_TYPE,`
 `WEIGHT_PATH,`
 `MAX_ITER_META,`
 `NUM_WORKERS,`
 `IMS_PER_BATCH,`
 `BASE_LR,`
 `BATCH_SIZE_PER_IMAGE,`
 `NUM_CLASSES)`

Parameters

`COCO_NAME: str`
 Image annotations exported to '.json' file in COCO format (<https://roboflow.com/formats/coco-json>)
`COCO_ANNOTS: str`
 Name of the annotated images folder
`MODEL_TYPE: str`
 Loading model zoo configuration, using ResNet and FPN (Feature Pyramid Networks) backbone
`WEIGHT_PATH: str`
 Loading pre-trained weights based on model zoo
 (https://github.com/facebookresearch/detectron2/blob/master/MODEL_ZOO.md) for instance segmentation
`MAX_ITER_META: int`
 Number of iterations
`NUM_WORKERS: int`
 Number of parallel data loading workers
`IMS_PER_BATCH: int`

Number of images per batch across all machines (depends on number of GPUs), each GPU will see 2 images per batch

BASE_LR: float
Hyperparameter that controls how much to change the model in response to the estimated error each time the model weights are updated (it has big impact for resulting model)

BATCH_SIZE_PER_IMAGE: int
Number of samples(images) that will be propagated through the network

NUM_CLASSES: int
Number of thing classes for R-CNN

4.1.3 Member Function Documentation

4.1.3.1 run() `def Modules.train.Trainer.run (self)`

Method for starting the training process including configuration for Detectron2

4.1.4 Member Data Documentation

4.1.4.1 BASE_LR `Modules.train.Trainer.BASE_LR`

4.1.4.2 BATCH_SIZE_PER_IMAGE `Modules.train.Trainer.BATCH_SIZE_PER_IMAGE`

4.1.4.3 COCO_ANNOTS `Modules.train.Trainer.COCO_ANNOTS`

4.1.4.4 COCO_NAME `Modules.train.Trainer.COCO_NAME`

4.1.4.5 IMS_PER_BATCH `Modules.train.Trainer.IMS_PER_BATCH`

4.1.4.6 MAX_ITER_META `Modules.train.Trainer.MAX_ITER_META`

4.1.4.7 MODEL_TYPE `Modules.train.Trainer.MODEL_TYPE`

4.1.4.8 NUM_CLASSES `Modules.train.Trainer.NUM_CLASSES`

4.1.4.9 NUM_WORKERS `Modules.train.Trainer.NUM_WORKERS`

4.1.4.10 WEIGHT_PATH `Modules.train.Trainer.WEIGHT_PATH`

The documentation for this class was generated from the following file:

- `train.py`

Index

- `__init__`
 - `Modules.train.Trainer`, 8
- `ap`
 - `run`, 6
- `arg`
 - `run`, 6
- `args`
 - `run`, 6
- `BASE_LR`
 - `Modules.train.Trainer`, 9
- `BATCH_SIZE_PER_IMAGE`
 - `Modules.train.Trainer`, 9
- `Canny`
 - `Modules.canny`, 3
- `COCO_ANNOTS`
 - `Modules.train.Trainer`, 9
- `COCO_NAME`
 - `Modules.train.Trainer`, 9
- `ext_modules`
 - `setup`, 7
- `help`
 - `run`, 6
- `IMS_PER_BATCH`
 - `Modules.train.Trainer`, 9
- `input_dir`
 - `run`, 6
- `main`, 2
 - `main`, 2
- `MAX_ITER_META`
 - `Modules.train.Trainer`, 9
- `MODEL_TYPE`
 - `Modules.train.Trainer`, 9
- `MODEL_USE1`
 - `Modules.config`, 5
- `MODEL_USE2`
 - `Modules.config`, 5
- `Modules.canny`, 2
 - `Canny`, 3
- `Modules.config`, 3
 - `MODEL_USE1`, 5
 - `MODEL_USE2`, 5
 - `TestMeta`, 3
 - `TestROI`, 4
 - `TrainMeta`, 4
 - `TrainROI`, 4
- `Modules.roi`, 5
 - `ROI`, 5
- `Modules.train`, 5
- `Modules.train.Trainer`, 7
 - `__init__`, 8
 - `BASE_LR`, 9
 - `BATCH_SIZE_PER_IMAGE`, 9
 - `COCO_ANNOTS`, 9
 - `COCO_NAME`, 9
 - `IMS_PER_BATCH`, 9
 - `MAX_ITER_META`, 9
 - `MODEL_TYPE`, 9
 - `NUM_CLASSES`, 10
 - `NUM_WORKERS`, 10
 - `run`, 9
 - `WEIGHT_PATH`, 10
- `name`
 - `setup`, 7
- `NUM_CLASSES`
 - `Modules.train.Trainer`, 10
- `NUM_WORKERS`
 - `Modules.train.Trainer`, 10
- `output_dir`
 - `run`, 6
- `required`
 - `run`, 7
- `ROI`
 - `Modules.roi`, 5
- `run`, 6
 - `ap`, 6
 - `arg`, 6
 - `args`, 6
 - `help`, 6
 - `input_dir`, 6
 - `Modules.train.Trainer`, 9
 - `output_dir`, 6
 - `required`, 7
 - `start_image`, 7
- `setup`, 7
 - `ext_modules`, 7
 - `name`, 7
- `start_image`
 - `run`, 7
- `TestMeta`
 - `Modules.config`, 3
- `TestROI`
 - `Modules.config`, 4
- `TrainMeta`
 - `Modules.config`, 4
- `TrainROI`
 - `Modules.config`, 4
- `WEIGHT_PATH`
 - `Modules.train.Trainer`, 10